

Moving and Handling Specials Index

Code	Product Title	Available Stock
	AA017288	Peacocks Snooooooze Sleep Small Wedge
	AA017635	Handicare Romedic Top Sheet ROMP 1515
	AA018090	Active Linx Easystand Evolve Standing Frame
	AA019160	Ross Return 7600 Patient Turner
	AA019569	Arjo Sara Plus
	AA019609	Handicare Turnsafe 2 Stand Aid
	AA020085	Mangar Leg Lifter and Airflow 12 Compressor
	SP291	Handicare Footstool

AA017288 Peacocks Snooooooze Sleep Small Wedge

- AA017288 - Auth Level:H



Peacocks Snooooooze Sleep Small Wedge

You must start an order to proceed.

Picture for illustrative purposes.

Product code: 9025

Wedge dimensions: 15x15x21 cm.

Fitting Instructions

This item cannot be fitted by GIS. GIS will deliver to the service users address for the assessor to fit to the new users requirements. Assessors must be competent in Moving & Handling Risk Assessment and safe Moving & Handling techniques. Care Plan instructions must be provided by the Assessor to user/ carers.

Select client to add 'Unsuitable' notes.

AA017635 Handicare Romedic Top Sheet ROMP 1515

- AA017635 - Auth Level:A



PRODUCT INFORMATION

Select client to add 'Unsuitable' notes.

Handicare Romedic Top Sheet ROMP 1515

You must start an order to proceed.

Product code: 1515

Fitting Instructions

This item cannot be fitted by GIS. GIS will deliver to the service users address for the assessor to fit to the new users requirements.

Assessors must be competent in Moving & Handling Risk Assessment and safe Moving & Handling techniques. Care Plan instructions must be provided by the Assessor to user/ carers.

If required, ask for support from your local Moving & Handling Keyworker. If you cannot identify a Keyworker, need training information or require specialist advice, contact Moving & Handling Specialists (MovingandHandling@glos.nhs.uk ; 0300 421 8175)

-AA018090 Active Linx Easystand Evolve Standing Frame

- AA018090 - Auth Level:A



PRODUCT INFORMATION

Select client to add 'Unsuitable' notes.

Active Linx Easystand Evolve Standing Frame

You must start an order to proceed.

Modular design makes the EasyStand Evolv the most versatile standing frame available

The base unit is functional for many users, however when more support is needed a wide variety of positioning options can be added to accommodate the most involved user

The Evolv enables full extension at the hips and knees for those with sufficient range of motion (ROM) and full weight-bearing on the long bones of the legs

For those with limited ROM, sit to stand provides an opportunity to stand in a partial weight-bearing position and potentially progress to full extension over time

Fitting Instructions

This item cannot be fitted by GIS.

GIS will deliver to the service users address for the assessor to fit to the new users requirements

-AA019160 Ross Return 7600 Patient Turner

- AA019160 - Auth Level:H



PRODUCT INFORMATION

Select client to add 'Unsuitable' notes.

Ross Return 7600 Patient Turner

You must start an order to proceed.

Model: ReTurn7600

Description Transfer platform, complete with wide rising ladder

Height 1210 mm
Total weight 26 kg
Baseplate, length 795 mm
Baseplate, width 672 mm
Baseplate, height 45 mm
SWL 205 kg

Fitting Instructions

This item cannot be fitted by GIS.

GIS will deliver to the service users address for the assessor to fit to the new users requirements.

-AA019569 Arjo Sara Plus

- AA019569 - Auth Level:H



Select client to add 'Unsuitable' notes.

Arjo Sara Plus

You must start an order to proceed.

Picture for illustrative purposes.

Model: SPL3030-F.

NOTE: This model has the scales option fitted.

Max safe working load kg 190 kg
Lift Weight (incl battery) kg 78,7 kg
Min total storage height mm 1195 mm
Min lifting height mm 906 mm
Max lifting height mm 1497 mm
External width legs closed mm 646 mm
External width legs open mm 1107 mm
Internal width leg closed mm 534 mm
Internal width leg open mm 1002 mm
Leg clearance (floor to top of legs) mm 116 mm
Chassis clearance (floor to bottom of chassis) mm 30 mm
Lift total length mm 1033 mm
Turning radius mm 1235 mm

Fitting Instructions

This item cannot be fitted by GIS. GIS will deliver to the service users address for the assessor to fit to the new users requirements.

Assessors must be competent in Moving & Handling Risk Assessment and safe Moving & Handling techniques. Care Plan instructions must be provided by the Assessor to user/ carers.

If required, ask for support from your local Moving & Handling Keyworker. If you cannot identify a Keyworker, need training information or require specialist advice, contact Moving & Handling Specialists (MovingandHandling@glos.nhs.uk ; 0300 421 8175)

-AA019609 Handicare Turnsafe 2 Stand Aid

- AA019609 - Auth Level:H



Select client to add 'Unsuitable' notes.

Handicare Turnsafe 2 Stand Aid

You must start an order to proceed.

Picture for illustrative purposes.

Name TurnSafe

Total weight 15 kg

Foot plate diameter 451 mm

SWL 150 kg

Unit 1

Spec_HeightOfLowerHandleUser - 800-1146 mm

Height of upper handle, for the caregiver - 912-1259 mm

Height of lower leg support - 185-531 mm

Lower leg supports adjustable in steps of - 50 mm

Fitting Instructions

This item cannot be fitted by GIS.

GIS will deliver to the service users address for the assessor to fit to the new users requirements.

-AA020085 Mangar Leg Lifter and Airflow 12 Compressor

- AA020085 - Auth Level:H



Select client to add 'Unsuitable' notes.

Mangar Leg Lifter and Airflow 12 Compressor

You must start an order to proceed.

Picture for illustrative purposes.

Lifting legs into bed can often prove difficult. Our inflatable Leglifter provides valuable assistance with lifting legs into bed using single handed care techniques. The Leglifter is attached to the bed underneath the mattress. It is straightforward to fit or remove and requires no permanent installation. When not in use the Leglifter sits discreetly against the side of the bed.

Model: MPCA140500.

Maximum User Weight: 135kg (21 stone).

Inflated Length: 124cm.

Inflated Width: 57cm.

Inflated Thickness: 10cm.

Deflated Thickness: 2 1/2cm.

Weight: 1.5kg.

Fitting Instructions

This item cannot be fitted by GIS.

GIS will deliver to the service users address for the assessor to fit to the new users requirements.

-SP291 Handicare Footstool

- SP291 - Auth Level:A



Select client to add 'Unsuitable' notes.

Handicare Footstool

You must start an order to proceed.

Picture is for illustrative purpose only

Product Code: 8090

Size: 320mm x 380mm x 100mm

SWL: 200 kg/440 lbs

Fitting Instructions

This item cannot be fitted by GIS. GIS will deliver to the service user's address for the assessor to fit to the new user's requirements.