



Walking Specials Index

Code	Product Title
AA014813	XXL Rehab Stand Tall
	
AA019677	GATE Bure StandTall Standard MANUAL Walker (Narrow)
	

- AA014813 XXL Rehab Stand Tall

- AA014813 - Auth Level:A



PRODUCT INFORMATION

Select client to add 'Unsuitable' notes.

XXL Rehab Stand Tall

You must start an order to proceed.

MAX USER WEIGHT: 325kg (37st)

DIMENSIONS:

- Adjustable Height: 1000-1300mm (39.3"-51.1")
- Width: 760mm (30")
- Internal Width: 710mm (28")
- Depth: 720mm (28.3")
- Front Wheel: 100mm (4")
- Rear Wheel: 120.5mm (4.8")

A multi-functional aid that can be used for independence to help users walk and mobilise from sitting to standing.

It has 2 built-in handgrips on the side which help the user to stand in small steps and in full control and brakes mounted to the bar which can be locked.

It has underarm-padding which is adjustable in width and angle.

Comes with IV holder and urine bag holder.

Fitting Instructions

This item cannot be fitted by GIS.

GIS will deliver to the service users address for the assessor to fit to the new users requirements.

- AA019677 GATE Bure StandTall Standard MANUAL Walker (Narrow)

- AA019677 - Auth Level:A



Manual

Select client to add 'Unsuitable' notes.

GATE Bure StandTall Standard MANUAL Walker (Narrow)

You must start an order to proceed.

Picture is for illustration purposes only

This version has 100mm caster wheels

Product Name - BURE S (Narrow) Manual

Product Code - RHGA-56-300

External Width - 66 cm

Length - 80 cm

Height - 91-133cm

Max User Weight 150 kg SWL

Fitting Instructions

This item cannot be fitted by GIS

GIS will deliver to the service users address for the assessor to fit to the new users requirements

Adjustments need to be carried out by a skilled health professional with the clinical knowledge to fit correctly for the user