


## Mattresses and Cushions Index

Code	Product Title
	<b>AA018483 Symmetrikit Foot Support Wedge</b>
	<b>AA018712 Invacare Sofform Bariatric Visco Mattress 193 x 152.5 x 15.2 cm</b>
	<b>AA019492 Lejrelet Oval Positioning Cushion</b>
	<b>AA019494 Lejrelet Positioning Tube 125 - Cushion</b>
	<b>AA020051 Dreama Single Bed Positioning System With Covers</b>
	<b>AA020053 Dreama Positioning Pad - Medium</b>
	<b>AA020054 Dreama Memory Foam Pad and Cover - Medium</b>
	<b>AA020300 LEJRELET Set With Tube 125 (5 Cushions)</b>

## -AA018483 Symmetrikit Foot Support Wedge

- AA018483 - Auth Level:A



### Symmetrikit Foot Support Wedge

You must start an order to proceed.

Note: picture is for illustration purposes only

Product Code: SA7660GB

Dimensions: Width: 70cm

Depth: 30cm

Height: 22cm

#### Fitting Instructions

**This item cannot be fitted by GIS**

**GIS will deliver to the service users address for the assessor to fit to the new users requirements**

Select client to add 'Unsuitable' notes.

## -AA018712 Invacare Soffform Bariatric Visco Mattress 193 x 152.5 x 15.2 cm

- AA018712 - Auth Level:E



### Invacare Soffform Bariatric Visco Mattress 193 x 152.5 x 15.2 cm

You must start an order to proceed.

Picture for illustrative purposes.

#### Fitting Instructions

This item cannot be fitted by GIS.

GIS will deliver to the service users address for the assessor to fit to the new users requirements

Select client to add 'Unsuitable' notes.

## -AA019492 Lejrelet Oval Postioning Cushion

- AA019492 - Auth Level:H



Select client to add 'Unsuitable' notes.

### Lejrelet Oval Postioning Cushion

You must start an order to proceed.

Picture for illustrative purposes.

Product code: 5000009.

Filling: The LEJRELET Oval is filled with 100% viscoelastic foam granules preventing decubitus.

Cover Polyester (Polyurethane coated).

Dimensions: 10 cm height x 25 cm width x 40 cm long.

Weight: 350g.

#### Fitting Instructions

**This item cannot be fitted by GIS.**

**GIS will deliver to the service users address for the assessor to fit to the new users requirements.**

## AA019494 Lejrelet Positioning Tube 125 - Cushion

- AA019494 - Auth Level:H



Select client to add 'Unsuitable' notes.

### Lejrelet Positioning Tube 125 - Cushion

You must start an order to proceed.

Picture for illustrative purposes.

Filling 100% viscoelastisk foam.  
Diameter 25 cm  
Length 125 cm  
Weight 2350 g  
Cover Polyester (Polyurethane coated)


#### Fitting Instructions

**This item cannot be fitted by GIS.**

**GIS will deliver to the service users address for the assessor to fit to the new users requirements.**

## -AA020051 Dreama Single Bed Positioning System With Covers

- AA020051 - Auth Level:A



Select client to add 'Unsuitable' notes.

### Dreama Single Bed Positioning System With Covers

You must start an order to proceed.

Picture for illustrative purposes.

Model: DC03.

User Weight Max: 110kg.

Double Cushion Size: LxWxH 900 x 95 x 150mm.

Double Cushion Size: LxWxH 900 x 195 x 150mm.

No. of Cushions Supplied: (single:double) 15:2.

Height of System (Inc. Base): 170mm.

Footprint: WxL 900 x 1900mm.

#### Fitting Instructions

This item cannot be fitted by GIS. GIS will deliver to the service users address for the assessor to fit to the new users requirements.

Set up of the sleep system needs to be carried out by a skilled health professional with the clinical knowledge to fit of the service user correctly.

## -AA020053 Dreama Positioning Pad - Medium

- AA020053 - Auth Level:A



Select client to add 'Unsuitable' notes.

### Dreama Positioning Pad - Medium

You must start an order to proceed.

Picture for illustrative purposes.

Model: Dreama Positioning Pad: DM09.

Made from comfortable and hygienic PU and available in 4 sizes. With or without memory foam pad for extra comfort.

#### Fitting Instructions

This item cannot be fitted by GIS. GIS will deliver to the service users address for the assessor to fit to the new users requirements.

## -AA020054 Dreama Memory Foam Pad and Cover - Medium

- AA020054 - Auth Level:A



### Dreama Memory Foam Pad and Cover - Medium

You must start an order to proceed.

Picture for illustrative purposes.

Model: DP02.

Memory foam pad for use with medium positioning pads DM09.

#### Fitting Instructions

This item cannot be fitted by GIS. GIS will deliver to the service users address for the assessor to fit to the new users requirements.

Select client to add 'Unsuitable' notes.



## -AA020300 LEJRELET Set With Tube 125 (5 Cushions)

- AA020300 - Auth Level:A



Select client to add 'Unsuitable' notes.

### LEJRELET Set With Tube 125 (5 Cushions)

You must start an order to proceed.

Picture for illustrative purposes.

The set includes:

LEJRELET Tube 125 is a 125 cm long cylindrical body support pillow. This support cushion can be used to support the body in the side-lying position, to calm, and provide comfort and tactile stimuli.

LEJRELET Oval is a small positioning pillow developed for micro-positioning. This small cushion can be used as an abduction cushion, lumbar support, or for relieving pressure on ears, elbows, hands, feet, and heels.

LEJRELET Pad High is a square and firm pressure relieving positioning cushion ideal for providing leg and knee support. Place one leg on top of the pillow to relax the hip and create a comfortable side-lying position.

LEJRELET Wedge is a slanted positioning cushion which can be used to create stable positions. The wedge-shaped cushion can be placed in front of the client's chest or at the back to create stability in the side-lying position.

LEJRELET Pad Low is a flat and moldable support cushion that can be rolled up or folded to suit the individual situation. For instance, fold the cushion and place it between the client's feet and the footboard of the bed, or roll it up and use it under the neck.

A practical bag to store and carry the set of positioning cushions.

#### Fitting Instructions

This item cannot be fitted by GIS. GIS will deliver to the service users address for the assessor to fit to the new users requirements.